

*Fits So Well...  
It Disappears*

AXIOM



The Axiom Ankle brace offers **staunch protection** to an injured ankle. But it **disappears** when the game begins.

*If an ankle brace is too flexible, you're not protected. If it is too stiff, it's distracting.*

- Stable... yet flexible
- Support... without restriction.
- Seamless fit in the shoe... and on the ankle.
- Low ankle/high ankle protection

***Prevent inversion & eversion.***

*You can FEEL the difference when you wear it — or can you?*

**B Bledsoe™**

*Engineered to Disappear. Engineered by Bledsoe.*



# 7 Steps to Make an Ankle Brace Disappear.



## 1. Find the Right Material

Uses a plastic alloy that prevents inversion & eversion, yet flexes with natural movement.

## 2. Fit, Don't Fight

Separate upper shells allow the material to naturally form to any ankle. The in-shoe foot plate easily heat-forms for a smoother fit.

## 3. Don't Forget the Syndesmosis

High ankle sprains are more complicated than low ankle sprains. The taller uprights support the high ankle. A second strap supports the syndesmosis.

## 4. Easy to Apply

An ankle brace must be quick. On and off is easy with innovative ski-buckles.

## 5. Remember Achilles' Defeat

Eliminate irritating pressure on the Achilles tendon, using a simple cut-out.

## 6. Eliminate Pressure Points

Center pivot point on the malleolus eliminates pressure and supports natural anatomical movement.

## 7. Easy to Customize

The unique footplate material is heat formable...at lower temperatures.

Faster. Easier. Safer.

*Enhance your  
DME program.*

Suggested Lcode: L1971

### ORDERING INFORMATION

	RIGHT	LEFT
SMALL	AS020203	AS020103
MEDIUM	AS020205	AS020105
LARGE	AS020207	AS020107

**Bledsoe™**

1.888.BLED SOE • LOCAL PHONE 972.647.0884 • INTERNATIONAL FAX 972.606.0649 • [www.BledsoeBrace.com](http://www.BledsoeBrace.com)



CP010150 Rev B 04/08

