

NEW!

Bledsoe EZ Set Hinged Boots

Hi-Top, Lo-Top, Mid-Calf, and PTB Styles



**Bledsoe
Brace
Systems™**



- **EZ Set Hinged Boot:** The Bledsoe EZ Set Hinged Boot provides complete range of motion settings from 20° dorsi-flexion to 40° plantar-flexion in 10° increments.
- **Enhanced mobility:** Designed to enhance mobility for sprains, strains, and various types of fractures.
- **Patient compliance assured:** Compliance straps provided to prevent patient tampering.
- **Reformable:** Aluminum boot shell and upright system that, unlike plastic boots, can be reformed to fit virtually any leg or ankle shape.
- **Comprehensive range of application:** The EZ Set Hinge may be applied to any premium Bledsoe Boot line; Hi-Top, Lo-Top, Mid-Calf or PTB (patella tendon bearing) Boots.
- **Optional Air Ankle / Heel Pad:** Allows for ability to adjust ankle pressure for increased comfort and provides pneumatic ankle support.

Superior support in a lightweight, EZ design



Hi-Top Boot style shown

Bledsoe EZ Set Hinged Boots

Restrictions on use

This device is offered for sale by or on the order of a physician or other qualified medical professional. This device is not intended for reuse on a second patient or for extended use beyond 180 days. This device is generally not intended for patients in excess of 250 lbs. (114 kg). This device is not intended to prevent injuries or reduce or eliminate the risk of injury.

Warning

Do not leave this device in the trunk of a car on a hot day or in any place where the temperatures may exceed 150°F (65°C) for any extended period of time. If any additional pain or symptoms occur while using this device, discontinue use and seek medical attention.



**Bledsoe
Brace
Systems™**

by Medical Technology, Inc.

2601 Pinewood Drive
Grand Prairie, TX 75051 USA

1-888-BLEDSOE
(1-888-253-3763)
1-800-527-3666

972-647-0884

Local Fax:
972-660-5495

International Fax:
972-606-0649

www.bledsoebrace.com

Indications

Hi-Top and Lo-Top: Indicated for post-injury or post-operative application to grade 1, 2, or 3 ankle sprains, and stable or internally fixed fractures of the foot or ankle. Also indicated for stable or internally fixed fractures of the middle to distal tibia or fibula when the optional tibial fracture bracing plate is used.

Mid-Calf: Indicated for post-injury or post-operative application to grade 1, 2, or 3 ankle sprains, and stable or internally fixed fractures of the foot or ankle.

PTB: Indicated for treatment of stable or internally fixed fractures or fractures with bridging callus of the middle to distal tibia or fibula when partial patellar tendon bearing is desired.

Contraindications

Hi-Top and Lo-Top: Contraindicated for unstable fractures and proximal tibial or fibular fractures.

Mid-Calf: Contraindicated for unstable fractures and proximal tibial or fibular fractures.

PTB: Contraindicated for unstable fractures and proximal tibial or fibular fractures.

Ordering Information EZ Set Hinged Boot

Note: Leg measurement taken from knee joint line to medial malleolus.

Shoe Size			Leg Measurement	Boot	Hi-Top Boot		Lo-Top Boot	
Men	Women	Euro	Length	Size				
1- 2.5	1-3.5	29-33	9"-11" (23 - 28 cm)	XS	w/cuff AL?72001	w/cuff AL?70001	w/bootie AL?72021	w/bootie AL?70021
2.5-5.5	3.5-6	33-37	11"-13" (28 - 33 cm)	S	w/cuff AL?72003	w/cuff AL?70003	w/bootie AL?72023	w/bootie AL?70023
5.5-9	6-9.5	37-41	13"-15" (33 - 38 cm)	M	w/cuff AL?72005	w/cuff AL?70005	w/bootie AL?72025	w/bootie AL?70025
9-12.5	9.5-13	41-46	15"-17" (38 - 43 cm)	L	w/cuff AL?72007	w/cuff AL?70007	w/bootie AL?72027	w/bootie AL?70027
12.5-17	13 +	46-49 +	17"-19" (43 - 48 cm)	XL	w/cuff AL?72009	w/cuff AL?70009	w/bootie AL?72029	w/bootie AL?70029

Note: When ordering replace the ? with one of the following: 0= ankle / heel pad 2= air ankle / heel pad

Note: Leg measurement taken from knee joint line to medial malleolus.

Shoe Size			Leg Measurement	Boot	PTB Boot		Mid-Calf Boot	
Men	Women	Euro	Length	Size	Left	Right		
1- 2.5	1-3.5	29-33	9"-11" (23 - 28 cm)	XS	w/cuff AL?73101	AL?73201	w/cuff AL?71001	
2.5-5.5	3.5-6	33-37	11"-13" (28 - 33 cm)	S	w/cuff AL?73103	AL?73203	w/cuff AL?71003	
5.5-9	6-9.5	37-41	13"-15" (33 - 38 cm)	M	w/cuff AL?73105	AL?73205	w/cuff AL?71005	
9-12.5	9.5-13	41-46	15"-17" (38 - 43 cm)	L	w/cuff AL?73107	AL?73207	w/cuff AL?71007	
12.5-17	13 +	46-49 +	17"-19" (43 - 48 cm)	XL	w/cuff AL?73109	AL?73209	w/cuff AL?71009	

Note: When ordering replace the ? with one of the following: 0= ankle / heel pad 2= air ankle / heel pad

Accessories:

- Rain, Cold Weather, or Night Cover
- Replacement cuffs or booties
- Tibial plates and Strap kits

Call for more information.

Bledsoe limited warranty

This device is warranted to be free from defects in material and workmanship for a period of 90 days. These warranties apply to devices that have not been modified or subjected to misuse, abuse or neglect.



For more detailed information on this brace including fitting and adjustment instructions, visit our Web site: www.bledsoebrace.com

For product information or questions pertaining to sales or service, please contact the sales representative in your area or Bledsoe Brace Systems.

