

Bledsoe Philippon K.A.F. Positioning Kit



**Bledsoe
Brace
Systems™**



The Bledsoe Philippon K.A.F. (Knee, Ankle, Foot) Positioning Kit was developed in collaboration with Dr. Marc Philippon, MD, Director of Sports Medicine/Hip Disorders at the University of Pittsburgh.

The Philippon K.A.F. Positioning Kit is designed for a dual function:

- **During the hip arthroscopic procedure:** The padded booties protect the skin, feet and ankles of the patient while pulling traction on the legs.
- **Following hip arthroscopy:** The booties, along with the knee band, hold the patient's feet and legs together to prevent hip rotation, adduction or abduction while allowing the patient to sleep comfortably.

The Philippon K.A.F. Positioning Kit is made from double-thick ultra-breathable polyurethane foam with a soft loop pile surface. The Philippon K.A.F. Positioning Kit pads and protects the patient in a comfortable manner

**Protects patients during
and after hip arthroscopy**

Bledsoe Philippon K.A.F. Positioning Kit

Restrictions on use

This device is offered for sale by or on the order of a physician or other qualified medical professional. This device is not intended for reuse on a second patient or for extended use beyond 180 days. This device is generally not intended for patients in excess of 250 lbs. (114 kg.) This device is not intended to prevent injuries or reduce or eliminate the risk of injury.

Warning

Do not leave this device in the trunk of a car on a hot day or in any place where the temperatures may exceed 150°F (65°C) for any extended period of time. If any additional pain or symptoms occur while using this device, discontinue use and seek medical attention.



**Bledsoe
Brace
Systems™**

by Medical Technology, Inc.

2601 Pinewood Drive
Grand Prairie, TX 75051 USA

1-888-BLEDSOE
(1-888-253-3763)

1-800-527-3666
972-647-0884

Local Fax:
972-660-5495

International Fax:
972-606-0649

www.bledsoebrace.com

Indications

This device is indicated to help limit hip rotation, adduction and abduction by positioning the knees, ankles and feet together following hip arthroscopy. It is further indicated to protect the patient's feet and ankles from high interoperative traction forces placed on the legs during hip arthroscopy.

Contraindications

None known.



Padded booties in the Philippon K.A.F. Positioning Kit protect the skin and nerves of the feet and ankles during surgery.

Knee and ankle bands in the Philippon K.A.F. Positioning Kit limit hip flexion following arthroscopic surgery.

Ordering Information Philippon K.A.F. Positioning Kit

Product	Size	Part Number
K.A.F. Positioning Kit	Universal	AT091005
Perineal Post Pad Kit	Universal	AT092000
Positioning with Post Pad Kit	Universal	AT092005

Bledsoe limited warranty

This device is warranted to be free from defects in material and workmanship for a period of 90 days for wraps, pads, straps and normal wear components. These warranties apply to devices that have not been modified or subjected to misuse, abuse or neglect.



For more detailed information on this brace including fitting and adjustment instructions, visit our Web site: www.bledsoebrace.com

For product information or questions pertaining to sales or service, please contact the sales representative in your area or Bledsoe Brace Systems.

