



**Bledsoe  
Brace  
Systems™**

## Bledsoe Mid-Calf Boot



Breathable foam wrap for maximum comfort

Two layers of pads contour to ankle shape, and provides total control and comfort

Strong aluminum alloy provides strength, yet allows formability to change widths

■ **Reformable:** Aluminum boot shell and upright system that, unlike plastic boots, can be reformed to fit virtually any leg or ankle shape.

■ **Low Profile:** Low profile rocker bottom helps eliminate pelvic tilt & alleviate low back pain.

■ **Patient comfort:** Ultra breathable foam wrap for maximum comfort.

■ **Fast, easy fitting:** Superior support without casting complications, and a shorter, lighter boot which requires less application time.

■ **Unique rocker design:** Allows patients to closely simulate a near normal gait without abnormal loads on any joints.

■ **EZ Set Hinged Boot:** EZ Set Hinged Mid-Calf boot is available that provides complete range of motion settings from 20° dorsiflexion to 40° plantarflexion in 10° increments.

■ **Optional Air Ankle / Heel Pad:** Allows for ability to adjust ankle pressure for increased comfort and provides pneumatic ankle support.

**Saves time  
because of  
its inherent  
simplicity**



EZ Set Hinge

# Bledsoe Mid-Calf Boot

## Restrictions on use

This device is offered for sale by or on the order of a physician or other qualified medical professional. This device is not intended for reuse on a second patient or for extended use beyond 180 days. This device is generally not intended for patients in excess of 250 lbs. (114 kg.) This device is not intended to prevent injuries or reduce or eliminate the risk of injury.

## Warning

Do not leave this device in the trunk of a car on a hot day or in any place where the temperatures may exceed 150°F (65°C) for any extended period of time. If any additional pain or symptoms occur while using this device, discontinue use and seek medical attention.



**Bledsoe  
Brace  
Systems™**

by Medical Technology, Inc.

2601 Pinewood Drive  
Grand Prairie, TX 75051 USA

1-888-BLEDSOE  
(1-888-253-3763)

1-800-527-3666

972-647-0884

Local Fax:  
972-660-5495

International Fax:  
972-606-0649

[www.bledsoebrace.com](http://www.bledsoebrace.com)

## Indications

Indicated for post-injury or post-operative application to grade 1, 2, or 3 ankle sprains, and stable or internally fixed fractures of the foot or ankle.

## Contraindications

Contraindicated for unstable fractures or for fractures of the proximal tibia or fibula.

## Bledsoe Mid-Calf Boot vs. other methods

Features	Bledsoe Boots	All other low profile walkers	All other profile walkers
Unique low profile rocker design allows for a near normal walking gait and alleviates lower back pain.	Yes	No (rocker too flat)	No (rocker curved, but too high)
Aluminum shell and uprights permit forming to fit most any size ankle or leg.	Yes	No	No
Ultra-breathable foam wraps optimize patient comfort and skin protection.	Yes	No	No
Open back design allows adjustment of foot and leg position over rocker to provide best walking gait.	Yes	No	No

## Ordering Information Mid-Calf Boot

Shoe size (preferred) or leg measurement taken from knee joint line to medial malleolus.

Shoe Size			Boot	Fixed Mid-Calf	EZ Set Mid-Calf
Men	Women	Euro	Size	Part Number	Part Number
1-2.5	1 - 3.5	29-33	XS	AL?55001	AL?71001
2.5-5.5	3.5 - 6	33-37	S	AL?55003	AL?71003
5.5-9	6 - 9.5	37-41	M	AL?55005	AL?71005
9-12.5	9.5 - 13	41-46	L	AL?55007	AL?71007
12.5-17	13+	46-49+	XL	AL?55009	AL?71009

**Note: When ordering replace the ? with one of the following: 0= ankle / heel pad 2= air ankle / heel pad**

### Accessories:

Rain, Cold Weather, or Night Cover - Replacement pads, cuffs, straps, plates, socks, and security ties  
**Call for more information.**

## Bledsoe limited warranty

This device is warranted to be free from defects in material and workmanship for a period of 90 days. These warranties apply to devices that have not been modified or subjected to misuse, abuse or neglect.



For more detailed information on this brace including fitting and adjustment instructions, visit our Web site: [www.bledsoebrace.com](http://www.bledsoebrace.com)

**For product information or questions pertaining to sales or service, please contact the sales representative in your area or Bledsoe Brace Systems.**

